

# Fume Hood Functionality Audit

## General Information

Lab / Location

Fume Hood ID / Number

Date of Audit

Auditor Name

## Inspection Checklist

Item	Pass	Fail	N/A
1. Sash is in good working condition	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Airflow monitor functioning	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. No obstructions inside fume hood	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Glass is clear and intact	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
5. Baffles are clean and unobstructed	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
6. Interior appears clean	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
7. Lighting is adequate	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
8. Electrical outlets working	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
9. Airflow within specification (if measured)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

## Comments / Actions Required

## Auditor Signature