Autoclave Sterilizer Maintenance Report Form

Facility/Department
Location
Autoclave Make & Model
Serial Number
Date of Maintenance
Performed By
Maintenance Checklist
Cleaned and inspected chamber
Cleaned and inspected chamber
Drain and refill water reservoir
Diam and Term water Teser von
⊫ Check door gasket and seals
Inspect electrical connections
Test cycle and verify operation
Lubricate moving parts
Calibrate temperature/pressure
Ohaa wati ana /Cammanta
Observations/Comments
Next Maintenance Due
Technician Signature
Technician Signature