Laboratory Water Bath Maintenance Checklist

General Information

Location	Model	
Serial No.	Date	
Checked by	Frequency	

Checklist

No.	Item	Yes	No	N/A	Comments
1	Water level adequate				
2	Temperature set correctly				
3	Temperature uniformity checked				
4	Cleanliness (interior/exterior)				
5	Drain and refill water				
6	Power cord/intact connection				
7	Functioning indicator lights/display				
8	Cover/lid damage check				
9	Leakage check				
10	Alarm function test (if applicable)				

Notes/Actions Taken:			
Technician Signature			
Date			
Supervisor Review			