

# Laboratory Fume Hood Inspection Checklist

Laboratory Name/ID

Fume Hood ID/Location

Date

Inspector

## Checklist

Item	Yes	No	N/A	Comments
Sash glass is intact and free of cracks	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<div></div>
Sash moves smoothly and stays in position	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<div></div>
Airflow monitor/alarm working	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<div></div>
Face velocity measured & within safe range	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<div></div>
Interior clean and free of clutter	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<div></div>
Lighting functioning	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<div></div>
Utility valves (gas, water) working properly	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<div></div>
Chemical containers capped and stored upright	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<div></div>
No obstructions to airflow at baffles/grill	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<div></div>
Certification current and visible	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<div></div>

Additional Comments

